



AUTHENTICATION SYSTEM BASED ON ECG

TEAM: SC_61



Team Members	
Name	ID
محمد صلاح احمد يوسف	
محمد محسن السيد عبدالرحمن محمد	
محمد رفعت ابو الحمد حمدان	
محمد سعيد محمد عبدالجليل	

Data Acquisition:

TB Diagnostic ECG Database

Data has been collected from <https://www.physionet.org/content/ptbdb/1.0.0/>

From 52 Healthy control subjects
we selected only 4 for our system

Data preparation:

First, we got the signals and removed isoline drift.

By using high pass filter and moving average filter to get the corrected filter.

(corrected=filtered - baseline)

Then do segmentation to extract ECG segments.

Then split data with ratio 80/20

And save one segment from each person in text file to be used for testing

Used Butterworth filter of order 4 and band (2,40)

Feature Extraction:

Then extract features using one methods:

1. Pan Tompkins algorithm for ECG points detection


Classification:

Then used 2 classifiers:


1. Support vector machines (SVM): parms=linear kernel
2. Linear Discriminant Analysis (LDA): parms= default

Classification Results:

	Classifier	
Feature extraction method	SVM	LDA
Fiducial features /Pan Tompkins	97.5 %	100.0 %
Fiducial features /11_points	95 %	97.5 %



SVM for 11 points



LDA for 11 points

	precision	recall	f1-score	support
sub_1	1.00	1.00	1.00	10
sub_2	0.91	1.00	0.95	10
sub_3	1.00	0.90	0.95	10
sub_4	1.00	1.00	1.00	10
accuracy			0.97	40
macro avg	0.98	0.97	0.97	40
weighted avg	0.98	0.97	0.97	40


SVM for

Pan


	precision	recall	f1-score	support
sub_1	1.00	1.00	1.00	10
sub_2	1.00	1.00	1.00	10
sub_3	1.00	1.00	1.00	10
sub_4	1.00	1.00	1.00	10
accuracy			1.00	40
macro avg	1.00	1.00	1.00	40
weighted avg	1.00	1.00	1.00	40

LDA

for Pan



SVM for Pan



LDA for Pan

Screen Shots:

Enter your Name[s]

Open file

Preprocessing

▾

Predict

Subject identified as sub_1 --> Access Allowed


Subject 1: p_156

Age: 17


Gender: Male


ECG date: 16/06/1992

signal



After preprocessing using fiducial_features





Enter your Name[s]

sub_5

Open file

sub_5

Preprocessing


optimized_Fiducial




Predict

Subject is unidentified --> Access
Denied

signal



After preprocessing using fiducial_features




Enter your Name[s]

Open file


Preprocessing ▾

Subject is unidentified --> Access Denied

signal



After preprocessing using pan_tompkins_11_poi



Architecture of the ECG identification

